




NONIN[®]

Disposable Sensors


Value Line of SpO₂ Sensors



NONIN's disposable sensors with PureLight[®] technology perform with patient comfort in mind.

- Ideal for extended monitoring or environments where risks of cross-contamination exist
- Micro-foam material for added comfort during continuous monitoring
- Repositionable tape with no adhesive transfer to patient
- Trusted PureLight technology, in a value-priced sensor

Sizes: Adult – Neonate



6000A - Adult
> 30 kg. / >66 lbs.

6000P - Pediatric
10 - 50 kg. / 22 - 110 lbs.

6000I - Infant
1 - 20 kg. / 2.2 - 44 lbs.

6000N - Neonate
<3 kg. / < 6.6 lbs.

PureLight Accuracy

PureLight[®] technology means accurate readings are assured.

Blood Oxygen Saturation
(%SpO₂ A_{rms}^{*}) 70 - 100% ± 2 digits for adults

Pulse Rate
±3% ±1 digit

www.nonin.com
800.356.8874
+1.763.553.9968



Please contact your authorized distributor for more information.

* A_{rms} encompasses 68% of the population.

© 2006 Nonin Medical, Inc.
Nonin[®], the Nonin Logo[®] and PureLight[®] are registered trademarks of Nonin Medical, Inc.
P/N 5583-000-01

NONIN[®] Reusable Soft Sensors

The Universal Sensor



For spot-checks or continuous monitoring, **NONIN's** soft sensors perform like two sensors in one, with patient comfort in mind.

- Quick and easy spot-check application
- Comfortable fit for continuous monitoring
- Durable design delivers consistent performance
- Suitable in many medical settings (EMS, Hospital, Sleep)
- Easy-to-clean with isopropyl alcohol
- Latex-free
- One-year warranty

Sizes: Adult – Neonate



8000SL - Large
Digit thickness 12.5 - 25.5 mm/
.5 and 1 inch



8000SM - Medium
Digit thickness 10 - 19 mm/
.4 and .75 inches



8000SS - Small
Digit thickness 7.5 - 12.5 mm/
.3 and .5 inches

PureLight Accuracy

PureLight[®] technology means accurate readings are assured.

Blood Oxygen Saturation
(%SpO₂ A_{rms}⁺) 70 - 100% ±2 digits for adults

Pulse Rate
±3% ±1 digit

www.nonin.com

800.356.8874

+1.763.553.9968

Please contact your authorized distributor for more information.

* A_{rms} encompasses 68% of the population. Not available for sale in Germany or Japan.

© 2006 Nonin Medical, Inc.
Nonin[®], the Nonin Logo[®] and PureLight[®] are registered trademarks of Nonin Medical, Inc.
P/N 5582-000-01

NONIN[®]

PureLight® Sensors – Simply Brilliant

Reusable SpO₂ Sensors

Nonin has designed an array of PureLight reusable sensors – available for patients of all sizes (adult, pediatric, & infant) and clinical presentations. The 8000 Series of reusable sensors are ideal for short-term monitoring, stress testing and spot checks – any situation where the risk of cross-contamination is low. A comfortable sensor, each is designed for a specific application and is form-fitted to decrease ambient light interference.



Reusable SpO₂ sensors:

- ✓ Proven PureLight LEDs for accuracy at critical SpO₂ levels
- ✓ Comfortable, self-aligning grips to assure proper placement
- ✓ Quick and easy application
- ✓ Easy-to-clean for convenience and reliability
- ✓ Variety of finger, earlobe and forehead sensors for various monitoring applications
- ✓ Cost-effective and durable, backed with our warranty

Sensor Type	1 meter	2 meter	3 meter	Extension Cable*
8000AP	✓		✓	✓
8000AA	✓	✓	✓	✓
8000K2	✓			✓
8000Q	✓			✓
8000R	✓			✓

*Optional Reusable Sensor Extension Cables are available in 1m, 3m, 6m and 9m lengths.

Reusable Flex SpO₂ Sensors with FlexiWrap®

Nonin's Flex Sensor System combines the convenience and performance of a disposable sensor with the economics of a reusable sensor. The reusable Flex Sensor, made of a durable silicone, is secured with a matching disposable FlexiWrap adhesive. This hybrid sensor system is ideal for numerous applications including extended monitoring and sleep studies. Available in sizes to fit neonates to adults, the Flex Sensor System is a great value for patient monitoring. Each sensor is packaged with 25 FlexiWraps (additional packages of 25 may be purchased separately).

Reusable Flex SpO₂ sensors with FlexiWrap:

- ✓ Proven PureLight LEDs for accuracy at critical SpO₂ levels
- ✓ Durability for repeated use
- ✓ Optimal performance in motion situations
- ✓ Cost-effective convenience
- ✓ Simple cleaning for ease-of-use
- ✓ Trim shape for comfort and contoured fit

Sensor Type	1 meter	2 meter	3 meter	Extension Cable*
8000J	✓		✓	✓
8008J	✓			✓
8001J	✓			✓

*Optional Reusable Sensor Extension Cables are available in 1m, 3m, 6m and 9m lengths.

Disposable SpO₂ Sensors

Nonin's Flexi-Form® II sensors are ideal for extended monitoring or situations where the risk of cross-contamination is high. The disposable Flexi-Form II sensors feature an inclusive cable, sensor and adhesive wrap to provide patients ultimate comfort, while providing a disposable solution at an attractive price. Each sensor has a simple placement illustration and an attached over-wrap that acts as a stabilizer for securing the sensor in position. These single-use sensors are available for patients of all sizes.

Disposable SpO₂ sensors:

- ✓ Proven PureLight LEDs for accuracy at critical SpO₂ levels
- ✓ Self-adhesive for extended monitoring and security
- ✓ Easy to apply for clinicians or patients
- ✓ Trim shape for comfort and contoured fit
- ✓ Can be cold cycle (EIO) sterilized
- ✓ Assortment pack (7000D) includes 4 adult, 2 pediatric, 2 infant, and 2 neonate

Sensor Type	1 meter	2 meter	3 meter	Extension Cable*
7000A	✓			✓
7000P	✓			✓
7000I	✓			✓
7000N	✓			✓

*Optional Reusable Sensor Extension Cables are available in 1m, 3m, 6m and 9m lengths.



7000A Adult Flexi-Form II
(> 30 kg / > 66 lb)



7000P Pediatric Flexi-Form II
(10 - 40 kg / 22 - 88 lbs)



7000I Infant Flexi-Form II
(2 - 20 kg / 4 - 44 lbs)



7000N Neonatal Flexi-Form II
(< 2 kg / < 4 lbs)

Adult Flex System



Infant Flex System



Neonate Flex System

